

Product Model:YM101/YM102/YM103

Version: 1.1  
Date: 2020-4-17

1 Product Introduction and Operation Guide

1.1 Front View




Figure 1 Front View of YM101/YM102/YM103

1.2 Operation Method

- Open the battery cover, and put the two AAA batteries into the battery compartment in correct polarities, then replace the cover;
- Press the bottom of the equipment and open the probe, then insert one finger into the probe;
- Press the button to turn the equipment on, and the measure interface will appear;
- After about 8 seconds, the measurement result can be read directly from the display screen;
- Before reading the parameters, make sure that stable numbers of the pulse oximeter interface has sustained more than 4 seconds;
- The equipment will turned off automatically within 8 seconds when the finger left the probe.

1

1.3 Battery Installation

- Put the two AAA batteries into battery compartment in correct polarities (Figure2).
- Push the battery cover horizontally along the arrow shown as right.




Figure 2 Battery Installation

1.4 Lanyard installation

- Pass the thinner end of the lanyard through the hanging hole;
- Pass the thicker end of the lanyard through the thinner end and tighten the lanyard (Figure3).

2



Figure 3 Lanyard Installation

1.5 Attention for Operation

- Before use check and confirm that the people or finger size were applicable;
- Before use check and confirm that the environment should be non-combustible material, as well as to avoid high or low temperature and humidity, but also need to pay attention to the following:
  - To avoid glare and direct sunlight exposure;
  - To avoid radiation infrared or ultraviolet radiation;
  - Avoid contact with the organic solvent, mist, dust, corrosive gases;
- The equipment should not be used at a location or limb tied with arterial canal or blood pressure cuff or receiving intravenous injection;
- The equipment may not work normally on microcirculation barrier patients, Warm or rub the finger, or re-position the equipment could improve the measurement.
- The ray between photo detector and light emitting diode should cross patient's arteriole.
- The patient should not use enamel or other makeup;
- Avoid to insert a wet finger into the probe.

3




Figure 4 Finger Placement Diagram

1.6 Functions and Menu Operation

1.6.1 The button operation rules

- Long-press functions include entering menus, activating item's submenu, confirming setting values, and exiting item's submenu; short-press functions are polling menu items and viewing the setting values of items. It should be noted that long-press means pressing the key

4

should be noted that long-press means pressing the key for about 2 seconds, and short-press means pressing the key for less than 0.5 second.

1.6.2 Menu Operation

Active the menu

- After the oximeter is turned on, long-press the power button to activate the menu, then short-press the button to view the setting values of each item. If the user wants to change the setting value of the item, long-press to enter the item's submenu, the parameter value starts to flash, short-press to traverse the parameter value until the parameter value required by the user is selected, long press to confirm and exit the submenu.

Item1. Setup the LED display brightness

- The first item is to setup the display brightness. Long-press the button to select a brightness level ranging from 1 to 3. The greater the value, the greater brightness of the display.




Item2.Setup the SpO2 Alarm Limits

- The second item is to setup the SpO2 alarm limits. For example: When SpO2 High limit is set to 96, an alarm will be issued when the spo2 value is higher than 92, and when the spo2

5

Spo2 low limit is set to 94, an alarm will be issued when the spo2 value is lower than 94.



Item3.Setup the PR Alarm Limits

- The third item is to setup the PR alarm limits. For example: When PR High limit is set to 130, an alarm will be issued when the PR value is higher than 130, and when PR low limit is set to 50, an alarm will be issued when the PR value is lower than 50.



Item4.Turn Alarm On/Off

- The forth item is long-press to turn Alarm on /off.



Item5.Check the software version

- The fifth item is to view the software version.

6

2 Specifications

2.1 Classification

- Type of protection against electric shock: II (Internally powered equipment)  
Degree of protection against electric shock: Type BF-Applied part  
Operating mode: Spot checking  
Degree of protection against hazards of explosion: IP22

2.2 Power Requirements

- Specification of battery: Two AAA (LR03)  
Operating current: 25-50mA

2.3 Physical Specifications

- Width\*Height\*Depth: 57×30×31 mm  
Weight: 28g (Bare machine)

2.4 Measurement Specifications

- SpO2 declared accuracy 70%~100%: ±2digits  
0% ~ 69%: unspecified  
30%~99%: 1%.

- PR declared accuracy : 25~250 bpm: ±3digits  
PR Resolution: 1 bpm

7

2.5 Environmental Specifications

- Temperature**  
Operating : +50~+104°F/+10~+40°C  
Storage/Transportation: -4~+104°F/-20~+60°C
- Humidity**  
Operating : 10~95%, noncondensing  
Storage/Transportation : 10~95%, noncondensing
- Atmosphere Pressure**  
Operating : 70~106kpa  
Storage/Transportation: 50~107.4kpa

2.6 Display

- Specification of battery: 1.5" LED Display;  
YM101: Red; YM102: Green  
YM103: White  
SpO2%, Pulse Rate, P1%, Bar Graph, Battery Indicator
- Display Type:  
Display Color:

Display content:

- Notes:  
1)The claim for oxygen saturation accuracy should be supported by clinical studies covering the entire claimed range.The fraction of inspired oxygen (FIO2) delivered to test subjects is varied to achieve a series of targeted steady-state saturation periods over the specified SpO2 accuracy range (e.g. 70 % to 100 %), then the SpO2 accuracy is calculated by comparing SpO2 readings of the pulse oximeter to the values of SaO2 determined with a Co-Oximeter.  
2)The clinical trial included 11 subjects, including 6 males and 5 females, with an age range of 18 to 46 years, the subjects skin color included dark black, medium black, light color and white.

8

3 Maintenance,Cleaning,Disinfection

3.1 Maintenance

- The equipment's design life expectancy is about 2 years, keep your equipment and accessories free of dust and dirt, and follow these rules:  
A. Please clean the equipment before use according to chapter 6.2;Remove the batteries inside the battery cassette if the equipment will not be operated for a long time;  
B. Replace the batteries in time when the battery voltage indicate lamps were empty;  
C. It is recommended that the equipment should be kept in a dry environment with no corrosive gases and good ventilation anytime. The moisture and high-light environments will affect its lifetime and even might damage the equipment.  
D. It is best to preserve the product in a place where the temperature is between -20 to 60°C and the relative humidity is less than 95%.  
E. The packed equipment can be transported by ordinary conveyance. The equipment not be transported mixed with toxic, harmful, corrosive materials.

WARNING

- No modification of this equipment is allowed.

9

3.2 Cleaning

- Your equipment should be cleaned on a regular basis. If there is heavy pollution or lots of dust and sand in your place, the equipment should be cleaned more frequently. Before cleaning the equipment, consult your hospital's regulations for cleaning the equipment. Recommended cleaning agents are:  
a) Mild soap (diluted).  
b) Ethanol (70%).  
To clean your equipment, follow these rules:  
a) Shut down the pulse oximeter;  
b) Clean the display screen using a soft, clean cloth dampened with a glass cleaner;  
c) Clean the exterior surface of the equipment and probe using a soft cloth dampened with the cleaner;  
d) Wipe off all the cleaning solution with a dry cloth after cleaning if necessary;  
e) Dry your equipment in a ventilated, cool place.  
To avoid damage to the equipment, follow these rules:  
CAUTIONS  
• Always dilute according to the manufacturer's instructions or use lowest possible concentration.  
• Do not immerse part of the equipment in the liquid.  
• Do not pour liquid onto the equipment or accessories.  
• Never use abrasive materials (such as steel wool or silver polish), or erosive cleaners (such as acetone or acetone-based cleaners).  
• If you spill liquid onto the equipment, contact us or your service personnel.

10

3.3 Disinfection

Clean the pulse oximeter before disinfecting it.The recommended disinfectant is ethanol 70%.Disinfection steps are the same as cleaning.

CAUTION

- Never use ETO or formaldehyde for disinfection.

3.4 Disposal

Dispose of the pulse oximeter in accordance with local environment and waste disposal laws and regulations.

4 Accessories

- One lanyard.
- Two AAA batteries(Optional).
- One user manual.
- One certificate card.

5 Troubleshooting

Trouble	Possible Reason	Solutions
The device can not be turned on	The batteries are drained away or almost drained away. The battery installation is incorrect. The device works abnormally.	Replace batteries. Install the battery over again. Please contact the product distributor

11

Trouble	Possible Reason	Solutions
The Spo2 and PR are not displayed normally	The finger size is too big or small Excessive ambient light User's blood perfusion is very low	Select the suitable size finger to measure Avoid the excessive ambient light irradiation Warm the finger and try again
The display is off suddenly	The device was set to shut down automatically in 8 seconds when there is no correct physiological signals The battery is almost drained away	Normal Replace batteries
The Spo2 and PR are not displayed stably	The finger is not inserted deep enough	Replace the finger and try again
	The finger is shaking or the body is moving	Try to keep still
	Not used in the work environment required by this manual	Please use in normal working environment
The device works abnormally.	The device works abnormally.	Please contact the product distributor

12

6 Appendix A EMC

The equipment complies with the requirement of standard EN 60601-1-2:2014 "Electromagnetic Compatibility – Medical Electrical Equipment".

1	Guidance and manufacturer's declaration – electromagnetic emission		
2	The model YM101/YM201 is intended for use in the electromagnetic environment specified below. The customer or the user of the model YM101/YM201 should assure that it is used in such an environment.		
3	Emissions test	Compliance	Electromagnetic environment – guidance
4	RF emissions CISPR 11	Group 1	The Model YM201 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
5	RF emissions CISPR 11	Class B	The Model YM201 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
6	Harmonic emissions IEC 61000-3-2	Not applicable	

13

Guidance and manufacturer's declaration – electromagnetic immunity			
The model YM101/YM201 is intended for use in the electromagnetic environment specified below. The customer or the user of the model YM101/YM201 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.

14

Electrostatic transient / burst IEC 61000-4-4	±2 kV for power supply lines 100 kHz repetition frequency ±1 kV for input/output lines	N/A	N/A
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode line-line	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % UT (100 % dip in UT) for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°	N/A	N/A
	0 % UT (100 % dip in UT) for 1 cycle at 0°	N/A	N/A

15

Guidance and manufacturer's declaration – electromagnetic immunity			
The model YM101/YM201 is intended for use in the electromagnetic environment specified below. The customer or the user of the model YM101/YM201 should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Power frequency magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 150 kHz to 90 MHz outside ISM bands	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the Models YM201, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter,

16

Radiated RF IEC 61000-4-3	10 V/m	10 V/m	Recommended separation distance $d \left[ \frac{3.5}{V_1} \right] \sqrt{P}$ $d \left[ \frac{3.5}{E_1} \right] \sqrt{P}$ 80MHz to 800MHz $d \left[ \frac{7}{E_1} \right] \sqrt{P}$ 800MHz to 2.7GHz
	80 MHz to 2.7 GHz		where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres(m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range <sup>a</sup>

17

Interference may occur in the vicinity of equipment marked with the following symbol: 			
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz to 10.15 MHz, 14 MHz to 14.2 MHz, 18.07 MHz to 18,17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.			

18

b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.			
c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the YM201 is used exceeds the applicable RF compliance level above, the YM201 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the YM101.			
d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

19

MANUFACTURER : Shenzhen Yimi Life Technology Co.,Ltd Add: 305, Building A, Tengbo Industrial Park, Changshangjiang Street, Longbei Village, Pingshan District, 518118, Shenzhen Tel: +86 755-86573112 Email: <a href="mailto:hnpzd@myspo2.com">hnpzd@myspo2.com</a> Web : <a href="http://www.myspo2.com">www.myspo2.com</a>			
EC REPRESENTATIVE Share Info Consultant Service LLC Repräsentanzbüro Add : Heerdter Lohweg 83, 40549 Düsseldorf E-mail : <a href="mailto:eu-rep@share-info.cn">eu-rep@share-info.cn</a>			

20

